

Nanoscale Disorder in Crystalline Materials

Stephan Rosenkranz^a, Stine Nyborg Ancona^a, Branton Campbell^b, Peter Chupas^c, and Ray Osborn^a

^a Materials Science Division, Argonne National Laboratory


^b Department of Physics and Astronomy, Brigham Young University

^c Advanced Photon Source, Argonne National Laboratory

Motivation

Many emerging phenomena of technological importance are strongly influenced by short range correlations on the 10-100Å length scale

Ferroelectrics




2005

Molecular Solids




Welberry 2003

Fast Ion Conductors




2005

CMR Oxides



Campbell 2001

Geometric Frustration




Lee 2002

Quasicrystals



Frey 2000

Metallic Alloys





Moss 2003


- Develop single crystal diffuse scattering to its full potential as one of the most powerful probes of both local distortions and short range correlations
- Use it in ongoing research program to study emerging complex phenomena

Neutron


The cross correlation method provides energy discrimination and efficient coverage of large volumes of reciprocal space.



- Simulations promise gains of up to a factor 100!




Reconstructed Scattering Function




Total Gain = Gain * #Channels
~ 0.3*300

Synchrotron X-ray



- Very efficient coverage of large volumes of reciprocal space
- Simplified corrections (absorption etc) when using high energy synchrotron X-rays (> 60 keV)



Rebinned data



- Example: diffuse scattering measured for nanoporous host-guest Prussian blue system



Cross Correlation

Prototype funded

SNS Letter of Intent

Prototype design, installation

Full proposal to SNS

Fully optimized instrument

~70

1985

2004

2005

2006

2007

Workshop

High-energy X-ray rotation method

Efficient rebinning, visualization

"Ultra-fast" detectors

Automatic refinement

S. Rosenkranz and R. Osborn, Neutron News 15, 21 (2004)